

广东环凯微生物科技有限公司

Guangdong Huankai Microbial Sci. & Tech. Co.,Ltd.

Add: No.788,Shenzhou Road,Guangzhou Science City,
Guangzhou,Guangdong Province,China

Tel: +86-20-32078333-8601 FAX:+86-20-32077703

E-mail: paulyanggz@yahoo.com.cn



Certificate of Analysis

Product HCM001
Plate Counts Agar 500g
Lot Number 230131A12
Expiration Date 2026.01.30

Appearance Characteristics	Specification	Results
Powder Colour	Pale yellow	Pale yellow
Dissolved Colour	Pale yellow	Pale yellow
Clarity	Clear	Clear

Physical Characteristics	Specification	Results
pH	7.0±0.2	7.17
Loss on drying	≤6%	4.26%

Biological Characteristics	Organism	Specification	Results
	<i>Escherichia coli</i> ATCC25922	PR≥0.7,Gray-White colonies	PR=0.83, Gray-White colonies
Recovery	<i>Staphylococcus aureus</i> ATCC6538	PR≥0.7, Yellow colonies	PR=0.89, Yellow colonies
	<i>Bcillus subtilis</i> ATCC6633	PR≥0.7,Gray-White colonies	PR=0.89,Gray-White colonies

Conclusion Qualified

Analyst 石小帆 Review 杨绪伟 Checker 