



GUANGDONG HUANKAI MICROBIAL SCI. & TECH. CO., LTD.

## Certificate of Analysis

HK/QR-Q/HKSJ 03-02

Product	Potato Dextrose Agar			N O. of Product	HCM050
Specs/ Type	500g	Lot	1110142	Expiration Date	2025.01.05
Date of Analysis	2022.01.10	Date of Report	2022.01.13	Standard Applied	Q/HKSJ03
Items		Specification		Results	
Appearance Characteristics	Powder Color		Pale Yellow powder	Pale Yellow powder	
	Dissolved Color,clarity		Pale Yellow, Translucent no precipitation	Pale Yellow, Translucent no precipitation	
Physical Characteristics	Moisture		≤6%	2,28%	
	pH		5.6±0.2	5.65	
Biological Characteristics	Function		Specification	Growth Results	
	Aspergillus niger ATCC16404		$P_R \geq 0.7$ , Mycelium white, black spores	$P_R=0.88$ , good growth mycelium white, black spores	
Conclusion	Qualified				
Analyst	姚华卓	Review	容艳芬	Checker	卢冠飞
				Seal	