

Guangdong Huankai Biotechnology Co. Ltd.
Report For Product Analysis

CHK/QR-JS-06-004

Product Name: Hydrogen Peroxide Test Strip (0-25mg/L)

Lot number: 3020171

Production Date: 2023.7.4

Report Date: 2023.7.4

Inspection Standard: Q/HK S03003

Analysis Result

Analysis Item	Standard Value	Analysis Result
Appearance	white	white
The Reaction Speed	complete in 15 seconds	qualified
Stability	The color is basically unchanged for 10 seconds after color development	qualified
Color Uniformity	color uniformity, no obvious dark spots	qualified
Color Rendering Accuracy	The difference between the measured value and the standard value should not exceed 1/2 color level	qualified
Test Blank	Test paper when testing 0mg/L standard point, Non-color	qualified
inspect the conclusion:		Qualified

Analyst: 顾源

Checker: 陈旭

Seal:

