



Guangdong Huankai Microbial Sci. & Tech.Co.,Ltd.

广东环凯微生物科技有限公司

Add: No.788,Shenzhou Road,Guangzhou Science City,
Guangzhou,Guangdong Province,China
Tel: +86-20-32078333-8654 FAX: +86-20-32078333-8803
E-mail: hugochiang@163.com

HCM117 C.L.E.D Medium

Usages:

For for the isolation ,enumeration and presumptive identification of microorganisms from urine.

Formulation (per liter):

Pancreatic digest of gelatin 4g

Pancreatic digest of casein 4g

Beef extract 3g

Lactose 10g

L-Cystine 0.128g

Bromthymol blue 0.02g

Agar 15g

Final pH: 7.3 ±0.2

Directions:

Suspend 36.1g in 1 L of distilled water , stirring heated to boiling until completely dissolved, dispensing flask, 121 °C autoclave for 15min.

Storage: Keep container tightly closed, store in a cool, dry place, away from bright light. Storage period of 3 years.

Specifications: 500g/bottle