

CERTIFICATE OF ANALYSIS



Version 1.0
Revision Date 2025-01-23
Date of First Issue 2025-01-23

GUANGDONG HUANKAI MICROBIAL SCI. & TECH. CO., LTD. REPORT FOR ANALYSIS OF DEHYDRATED MEDIA

Inspection Number: 250233

HK/QR-Q/HKSJ 03-02

Product Name	Tryptic Soy Broth		Product Code	HCM057		Date of Production	2025.02.20	
Specs/ Type	500g		Lot	250220A12		Valid Till	2028.02.19	
Date of Analysis	2025.02.21		Date of Report	2025.02.22		Standard Applied	Q/HKSJ 03	
Test Parameters			Specification			Results		
Appearance Characteristics	Powder Color		Light yellow			Qualified		
	Dissolved Color		Light yellow			Qualified		
	Clarity		Clear without precipitation			Qualified		
Physical Characteristics	pH		7.3±0.2			7.42		
	Moisture		≤6.00%			2.58%		
Biological Characteristics	Function	Quality Control Strain Name	Specification		Growth Results			
			Inoculum (cfu/dish)	Growth Characteristics				
		<i>Bacillus subtilis</i> ATCC6633	<100	Broth is cloudy	Characteristic conformed			
		<i>Candida albicans</i> ATCC10231	<100	Broth is cloudy	Characteristic conformed			
	<i>Aspergillus niger</i> ATCC16404	<100	Broth is cloudy	Characteristic conformed				
Conclusion	Qualified							
Analyst	容乾芬		Reviewer	杨绪伟		Approver		