



Guangdong Huankai Microbial Sci. & Tech.Co.,Ltd.

广东环凯微生物科技有限公司

Add: No.788,Shenzhen Road,Guangzhou Science City,
Guangzhou,Guangdong Province,China

Tel: +86-20-32078333-8654 Website:www.hkmbio.com

E-mail: hugochiang@163.com

022039P1 Mac Conkey Broth(Granular Culture Media)

Usages:

For cultivation of *Escherichia coli*.

Principle:

Peptone provide nitrogen, vitamins and growth factors; sodium chloride to maintain osmotic equilibrium; lactose into fermentable sugars.

Formulation (per liter) :

Pancreatic Digest of Gelatin	20 g
Lactose Monohydrate	10g
Dehydrated Ox Bile	5g
Bromocresol Purple	10mg
final pH	7.3 ±0.2

How to use:

1.Suspend 35g in 1 L of distilled water , stirring heated to boiling until completely dissolved, dispensing flask, 121 °C autoclave for 15min.

2.Diluted and treated samples.

Storage: Keep container tightly closed, store in a cool, dry place, away from bright light. Storage period of 3 years.

Specifications: 250g/bottle