



Guangdong Huankai Microbial Sci. & Tech.Co.,Ltd.

广东环凯微生物科技有限公司

Add: No.788,Shenzhou Road,Guangzhou Science City,
Guangzhou,Guangdong Province,China
Tel: +86-20-32078333-8654 FAX: +86-20-32078333-8803
E-mail: hugochiang@163.com

HCM049 Sabouraud Dextrose Agar(SDA)

Usages:

For the cultivation of *Candida albicans*.

Principle:

Glucose as carbon source; peptone provide nitrogen; agar as medium coagulant

Formulation (per liter) :

Glucose 40g

Peptone 10g

Agar 15g

Final pH 5.6 ±0.2

How to use:

1.Suspend 65g in 1L of distilled water , stirring heated to boiling to completely dissolve ,autoclave at 121 °C for 15 minutes.

2.Diluted and treated samples.

Storage: Keep container tightly closed, store in a cool, dry place, away from bright light. Storage period of 3 years.

Specifications: 500g/bottle

Reference:

1. "People's Republic of China Pharmacopoeia" 2010 edition, Appendix XI J microbial limits.

2 ISO 18416-2007 Cosmetics ---- Microbiology ----- Detection of *Candida albicans*

Test national standard pharmaceutical industry cleanroom 3.GB/T 16294-2010 People's Republic of China (area) Settlement Act