

Product Manual

【Product Name】

Novobiocin (Additive for Modified EC Novobiocin Broth)

【Product Code and Packing Specification】

Code	Type	Specifications
SR0030	Lyophilized reagents	10 vials/box

【Product Usage】

Each vial is added to 225 ml (028112) to prepare modified EC neomycin enrichment broth.

【Formulation】

Formulation (per vial)	Content
Novobiocin	4.5mg

【Instructions】

1. How to open the vial: There is an arrow mark on the aluminum cap of the vial. Before opening the vial, disinfect the surface of the vial with 75% alcohol cotton. Under sterile conditions, open the aluminum cap in the direction of the arrow on the aluminum cap, tear off the aluminum cap, and open the vial rubber stopper.
2. How to use the reagent: Add 1mL of sterile saline to each vial to completely dissolve the reagent, and then add each vial to 225mL (028112) to prepare the modified EC neomycin enrichment broth.

【Storage Conditions and Shelf Life】

Store at 2-8°C, valid for two years.

【Waste Disposal】

After testing, the contaminated items need to be autoclaved at 121°C for 30 minutes before processing.

【Executive Standard】

Q/HK 0709 Microbiological Detection Kit

【Explanatory Version】

May 14, 2019

Guangdong Huankai Microbial Sci. & Tech. Co.,Ltd.

Website: www.huankai.com

Address: No.788,Shenzhou Road,Guangzhou Science City,
Guangzhou,Guangdong Province,China

Post code: 510663

Fax: 020-32078333-8619

Sales hotline: 020-32078333 ext. 8602 (extension)

Technical hotline: 020-32078333 ext. 8877, 8876 (extension)